

SMAS FOCUS V

/ Skin Tightening · Face Lifting · Body Slimming · Vaginal Tightening /




5D HIFU



One Machine, Multi-Functional

Only need to replace the handle, you can upgrade the product function at any time.

How it works




HIFU-powered transducers deliver concentrated waves of ultrasound into multiple layers of skin and body tissues without making direct contact with surrounding areas.

Transferring stabilized temperatures of 65~75°C into the targeted zones, the SMASFOCUS V aims to induce a process of neocollagenesis in the face while tightening body tissues for contouring purposes


FOR LIFTING & TIGHTENING




1 The dermis layer of the skin naturally synthesizes less collagen over the course of time.



2 Focused ultrasound waves deliver energy directly within the skin layer causing a shrinking effect.



3 Stimulated by the focused energy, skin collagen starts to regenerate.




4 Through the regeneration of collagen, wrinkles are eliminated in the process of both firming and tightening the skin.

FOR SLIMMING & CONTOURING



1 Fat cells and cellulite is formed within the body.



2 The SMASFOCUS V energy destroys the fat cells that have been targeted.




3 As a natural healing response of the body the cells regenerate joining together to fill the gap.



4 In result, the layer of fat naturally exits the body delivering a smooth, contoured body line.

FOR VAGINAL TIGHTENING



Product Introduction





Vmax single dot handle



1.5mm, 3.0mm, 4.5mm

Eyebrow, Wrinkles around eyes,
Wrinkles around mouth



2D hifu handle



1.5mm, 3.0mm, 4.5mm

Check, Forehead, Double chin,
Neck wrinkles, V-line forming



3D hifu handle



1.5mm, 3.0mm, 4.5mm

Check, Forehead, Double chin,
Neck wrinkles, V-line forming



6.0mm, 8.0mm, 10.0mm, 13.0mm, 16.0mm

Abdomen, Inner thighs, Outer
thighs, Hips & Buttocks, Inner
arms, Breast



Liposonix handle



8.0mm, 13mm

Body slimming and fat reduction




Vaginal handle



3.0mm, 4.5mm

Vaginal tightening



Easy to upgrade


Each machine is pre-loaded with a software system with all handle functions before shipment. When users need to upgrade the machine functions, they only need to purchase the handle with that function, plug the handle into the machine, and the software system can automatically identify the new handle to realize machine function upgrade.




2D HIFU



3D HIFU



Vmax



Liposonix



Vaginal

Summary of Advantages

L Procedure Time
30 - 45 Minutes

♥ Full recovery
Immediate

+ Number of Treatments
1 every 12 - 18 months

✓ Back To Work
Immediate

📎 Software Perloaded
All functions' software perloaded already

⚠ Risks & Complications
Possible mild redness, swelling and tenderness post treatment

🪄 Multifunctional
5 handles optional

📈 Upgrade
By inserting new handles only

i Duration Of Results
12 - 18 months



Before and after



Before

After



Before

After



Before

After



Before

After



Before

After



Before

After

Specifications

Power 1200 W

Frequency 4 ~ 7MHz

Net Weight 25kg

Dimension 400*350*1060mm

Voltage AC 100-240v 50/60Hz

