




Spheresculpt therapy machine



- Body Contouring
- Lymphatic Drainage
- Skin Tightening
- Pain Relief
- Vascular Action
- Muscle Toning



Principles of treatment



Spheresculpt Therapy produces a vibration and pressure combination that performs in effect gives the skin a 'workout'. This generates the drainage of fluids, the re-compacting of skin tissues, the removal of the "orange peel" effect from under the skin's surface. It also helps microcirculation which can help reduce inflammation and improve in muscle tone.

It helps improve vascularisation which in turn supports the production of collagen and elastin. It increases oxygen delivery to help nourish and brighten tissue from within. It tones the muscles helping reduce eliminate the appearance of expression wrinkles, combat tissue sagging, and generally lift the complexion and facial structure.

Range of treatment



Double Chin



Axilla



Inner Thighs



Flanks




Bra Line



Abdomen



Love Handles



Below Buttocks



Outer Thighs

Body treatment

/ Areas of application /

Cellulite

Contracted Muscles

Swollen legs

Buttocks

Abdominal fat

Post-pregnancy treatment

Orange Peel skin

Oedema

Arms

Ankles

Breast tightening

- REDUCES LYMPHATIC FLUID RETENTION
- REDUCES ORANGE PEEL SKIN
- IMPROVES VENOUS CIRCULATION
- IMPROVES BODY SHAPE
- TONES MUSCLES



FACIAL TREATMENT

/ Areas of application /

Mouth

Forehead

Cheekbone

Neck and lower neck area

Eyes

Cheeks

Glabella

BreastPre and post Injectable treatments



- SMOOTHES WRINKLES
- IMPROVES VASCULARIZATION AND OXIGENATION
- LIFTS AND FIRMS THE SKIN
- RELAXES CONTRACTED MUSCLES (EXPRESSION LINES)
- STIMULATES COLLAGEN PRODUCTION



Advantages

Spheres

The maniple roller consists of 11 rows with 5 balls in each, 55 balls in total. The spheres are made of hypoallergenic silicone, which is approved for use in the food and medical industry

Fast & Effective

Visible results after first visit. Skin quality improvement. Fat reduction up to 2 clothing sizes

Three kinds of spheres

Three kinds of spheres of different materials can be easily applied to different treatment items



Noise Level

With noise level less than 49 dB working on the maximum speed, treatment is more comfortable for both, client and specialist.

Maniple Design

The treatment is Smooth and soft impact and leaves no bruises and leaves no bruises. No special clothing is needed.

Usability

With ideal weight of 2,350 g there is no need to apply additional pressure to the maniple during the treatment.

Before and After



SPECIFICATIONS

Screen	15.4 inch
Rotating speed	450rpm
Output frequency	40-300Hz
Output power	300w
Voltage	AC110V/220V
Fuse	10A
Dimension	420x440x1220mm
Net weight	40kg



MED - LUX
2024