


# MED-LUX

Accentuate your natural beauty and renew your confidence with Vshape



# LOVE YOUR SHAPE

## SAY GOODBYE TO CELLULITE




Whole body treatment

Safe and Reliable

Effective and Comfortable






Infrared laser directly heats the targeted area

A vacuum elevates the target tissue, bring it closer to the energy source

Bi-polar RF heat from epidermis to subcutaneous at 15mm depth, can deliver thermal energy continuously to adipose tissue and dermal collagen fibers.


Roller massage accelerates blood cycle and lymphatic drainage.



## Combines five different technologies


The V-shape combines five different technologies including IR (Infrared laser), Bi-polar RF, massage roller, vacuum and ultrasound cavitation. The combination of IR, vacuum and RF cause deep heating of the fat cells, their surrounding connective tissue and the underlying dermal collagen fibers. This type of efficient heating and vacuum stimulates the growth of new and better collagen and elastin which results in localized reduction in skin laxity, body volume, and an overall improvement in skin structure and texture. The massage roller causes an immediate increase in circulation and lymphatic drainage, both essential components for healthy skin structure.






## Face and Submental

- Skin rejuvenation & skin tightening (Laxity, wrinkles, fine lines)
- Face & neck contouring (submental shaping)
- Skin resurfacing




## Arms (bingo wings) & Back

- Circumferential reduction
- Skin tightening (Laxity)
- Stretch marks



## Stomach (abdomen & love handles)


- Circumferential reduction
- Body contouring
- Skin tightening (Laxity)
- Stretch marks (red and white)



## Legs (thighs & buttocks)

- Circumferential reduction
- Cellulitis treatment
- Body contouring
- Skin tightening





Front

Back





40mmx66mm

### INFRARED VACUUM ROLLER RF HANDLE (BIG)

Used for cellulite removal  
skin tightening on the abdomen, buttocks, waist, thighs and back.

30 seconds of treatment can  
increase the temperature of the  
treatment area to 43.5°C





**Front**

**Back**

Clean And Hygienic


Applicator head replaceable design



30mmx44mm

**INFRARED VACUUM ROLLER RF HANDLE (SMALL)**

Used for cellulite removal and skin tightening on arms, calves and back.



Front

Back



16mmx25mm




9mmx13mm




4.5mmx9mm

## FACIAL VACUUM RF HANDLE


Equipped with three different sizes of treatment heads, used for facial and neck skin tightening and lifting, eye wrinkles removal, eye bags and dark circles treatment.





Front

Back




The handles adopt 3.5-inch touch screen, and the treatment parameters can be adjusted by the handle during treatment, which is convenient for operation



## VACUUM ULTRASOUND CAVITATION HANDLE


Combines vacuum and ultrasound, for local fat reduction in the abdomen, buttocks, thighs and other areas.

## AVERAGE REDUCTION OF ABDOMEN WIDTH



Based on numerous clinical medical treatment statistics, it shows 4-6 cm waist width reduction after one cycle of treatment, and this reduction keeps going on for 2-3 future months.

## BEFORE AND AFTER








### Japan SMC Solenoid Valve

High speed, long life, no noise



### Korea DKC Motor

Stable output, long life



### Italy Imported Air Pump

Large flow, low noise



### America Coherent IR Laser Diode

High power, long life

Model Number	MED-LUX
Power of RF	800W
Frequency of RF	2Mhz
Infrared laser Power	20w
Infrared laser wavelength	940nm
Ultrasound power	200W
Vacuum Indensity	10-90Kpa
Dimensions	532mm x 537mm x 1130mm
Weight	50kg

